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thermo scientific

 **iNORGANIC™**  
VENTURES

 **APPLIED SPECTRA**  
*Transforming the way the world does chemistry*

 **AL**  
**CHEM-LAB**

redoks  
LAB

**ThermoFisher**  
SCIENTIFIC

# ThermoFisher SCIENTIFIC

Type	Public
Traded as	NYSE: TMO <a href="#">↗</a> S&P 100 component S&P 500 component
Industry	Laboratory equipment, biotechnology, pharmaceutical, healthcare
Founded	1956 <sup>[1]</sup>
Headquarters	Waltham, Massachusetts, U.S.
Area served	Worldwide
Key people	Marc N. Casper (President & CEO) Mark Stevenson (COO) Ryan Snyder (CIO)
Products	Analytical/other equipment and instruments, laboratory reagents and consumables, science software and services —for research, discovery, analysis, and manufacturing (incl. pharmaceutical and diagnostic products)
Revenue	▲ US\$39.21 billion (2021) <sup>[1]</sup>
Operating income	▲ US\$10.03 billion (2021) <sup>[1]</sup>
Net income	▲ US\$7.73 billion (2021) <sup>[1]</sup>
Total assets	▲ US\$95.12 billion (2021) <sup>[1]</sup>
Total equity	▲ US\$40.79 billion (2021) <sup>[1]</sup>
Number of employees	c. 130,000 (Dec 2021) <sup>[1]</sup>

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## Sample Requirements Criteria

Criteria	Flame AA	GFAA	ICP-OES	ICP-MS
<b>Measurement Range</b>				
high > 10%			X	
1 - 10 %	X		X	
ppm	X		X	X
high ppb	X	X	X	X
low ppb		X	X	X
ppt		X		X
<b>Number of samples</b>				
Few	X	X		
Several	X		X	X
Many			X	X
<b>No Elements per Sample</b>				
Single	X	X	X	X
Few (2-5)	X		X	X
Intermediate (5-10)			X	X
Many			X	X
<b>Sample Matrix</b>				
< 3%	X	X	X	X
3-10 %	X	X	X	
> 10%		X	X	

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thermo**scientific**

APPLICATION NOTE 43121

## The analysis of gold using the Thermo Scientific ICE 3000 Series Atomic Absorption Spectrometers

### Authors

Dr. Jianfeng Cui,  
Thermo Scientific,  
Bremen, Germany

### Keywords

Cyanide leach, Fire Assay,  
Flame Atomic Absorption, Gold

### Introduction

Gold is a brilliant and beautiful metal, often associated with wealth, power and luxury. It is malleable and ductile, highly resistant to corrosion and tarnishing and sought after around the globe. Its uses range from jewellery and decoration, to storage of monetary value; to dentistry. It can also be found in food, electronics, and medicine.

Gold mining is the process of removing gold from the ground or streams, followed by multiple processing and extraction stages to obtain the final desired product. Due to its importance and high value, gold prospectors and metallurgists need to accurately determine the concentration of gold in mined samples. As geological samples are typically not homogenous, drying, grinding and pulverizing of samples is required before dissolution and analysis. Analysts require rapid and robust analytical methods on instrumentation capable of handling high levels of dissolved solids, yet providing high sensitivity, accuracy and stability. Flame atomic absorption spectrometry is the ideal technique to meet these analytical needs. This application note presents two methods for analysing gold. The first uses an aqua regia sample digestion method, while the second follows a cyanide leach protocol. Quality Control checks are used to ensure the long term stability of each method.

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#### Author

Dhinesh Asogan,  
Application Specialist,  
Thermo Fisher Scientific,  
Bremen, Germany

#### Keywords

Catalyst, Electronics, Jewelry,  
Metallurgy, Microwave digestion,  
PGE, Platinum group elements

#### Goal

To demonstrate the ability of the Thermo Scientific iCAP PRO XP ICP-OES to perform analysis of Platinum Group Element (PGE) analytes in high concentration metal matrices with excellent accuracy, where traditionally a radial instrument would be used.

#### Introduction

The platinum group elements, or PGEs as they are sometimes referred to, are a collective of metals grouped together in the middle of the periodic table. The six elements normally included in the group are all transition metals with some similar chemical and physical properties. Generally, this group is assumed to include platinum, palladium, rhodium, iridium, ruthenium and osmium (see Figure 1); all have become valuable and extremely useful to industry (and civilization in general) because of their elemental properties, rarity and beauty.

A wide range of PGE alloy compositions are used in electronic applications such as low voltage and low-energy contacts, thermocouples, furnace components and electrodes. The chemical and petrochemical industries use the PGEs' extraordinary catalytic properties in a wide range of processes including refining crude oil to produce many synthetic organic chemicals.<sup>1</sup>

SelectScience® 24/7  
The Fastest Way to Expert Opinion™

Search products, reviews, methods, videos...

News & Advice | Products & Reviews | Write a Review | Special Features | The Scientists' Channel | Login | Join Free

General Laboratory | Life Sciences | Pharma | Clinical | Analytical Chemistry | Applied Sciences | Automation & Informatics


Products / iCAP PRO Series ICP-OES

## iCAP PRO Series ICP-OES by Thermo Fisher Scientific

Manufacturer Thermo Fisher Scientific | Model: IQLAAGGAARFARDMZZZ | Available Worldwide

★★★★★ 4.9 / 5.0 | 7 reviews | Write your own review

Analyze the most challenging samples

iCAP PRO Series ICP-OES

Great instrument!

Yu Dong  
lubrizol

### Why request a quote through SelectScience?

- ✓ **Easy**  
It's quick and simple to do.
- ✓ **Fast**  
Your inquiry will be delivered straight to the manufacturer.
- ✓ **Free**  
You're under no obligation.
- ✓ **Secure**  
We only pass your details on to trusted suppliers at your request.
- ✓ **Save time**  
Submit your details once and make multiple inquiries.

**Request Pricing**

Turn On Instant Requests

**ThermoFisher SCIENTIFIC**

Receive your quote directly from the manufacturer.

**Quick Links**

- Request Product Information
- Contact Supplier

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### Authors

Daniel Kutscher,  
Simon Lofthouse,  
Simon Nelms and  
Shona McSheehy Ducos

### Keywords

Isobaric interferences,  
Kinetic energy discrimination,  
Reaction gas,  
Single quadrupole ICP-MS,  
Triple quadrupole ICP-MS

### Goal

To show the complete removal of isobaric interferences such as  $^{87}\text{Rb}$  on  $^{87}\text{Sr}$ , or  $^{204}\text{Hg}$  on  $^{204}\text{Pb}$  utilizing the Thermo Scientific ICAP TQ ICP-MS.

### Introduction

Unresolved spectral interferences may lead to biased results in ICP-MS. Single quadrupole ICP-MS typically uses a comprehensive interference removal mechanism with an inert collision gas (helium) and kinetic energy discrimination (KED). This approach, is based on the difference in size between an analyte and polyatomic interference and as such is effective for removal of many common interferences. In some cases though, two elements may share isotopes with identical mass number. These isotope overlaps are commonly referred to in ICP-MS as isobaric interferences. As the elemental ions are of a similar size, KED is not able to resolve these interferences. In addition, as the difference in mass between overlapping isotopes is extremely low, they cannot be spatially resolved using either quadrupole based or high resolution ICP-MS. Often, correction of isobaric interferences is accomplished using mathematical methods, however this procedure may lead to increased measurement uncertainties. Isobaric interferences may be resolved using reactive gases (e.g.  $\text{O}_2$  or  $\text{NH}_3$ ), in the cases where one element forms a different product ion to the other. However, side reactions may create new unwanted interferences with other ions extracted from a sample, so that full interference removal may not be possible.

### **Laser Induced Breakdown Spectroscopy (LIBS)**

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### Chemical Imaging of Ruby-in-Zoisite Gem Rock Using J200 LIBS Instrument

Elemental mapping for ruby-in-zoisite gem rock using Applied Spectra's J200 LIBS Instrument demonstrates how LIBS, with the ability to measure lighter and organic elements, is a highly complementary elemental imaging technique for XRF/XRD, electron microprobe, and LA-ICP-MS.


### Geochemical Fingerprinting of Coltan Minerals Using LIBS

LIBS spectra contains rich composition information that can help analysts identify the origin or source of geological samples. Read this application note to discover how Applied Spectra's J200 LIBS Instrument identifies coltan mineral ores important for electronic manufacturing from conflict zones.

### Rapid Lead (Pb) Analysis of Thin Solder Plating on Semiconductor Leadframes Using LIBS

Read this application note to understand how powerful the J200 LIBS Instrument is for sensitive micro-analysis of electronic components and how LIBS compares with traditional elemental analysis techniques such as XRF and ICP-OES.

### Innovative Elemental Mapping of Geological Minerals with Applied Spectra's J200 Tandem LA-LIBS

To understand the capabilities of the J200 Tandem LA - LIBS instrument in combination with ICP-MS, a rare earth element (REE)-rich mineral was analyzed and the elemental composition over a 16 mm<sup>2</sup> area was mapped using contour plots.

### Rapid Analysis of Mining Samples Using Laser Induced Breakdown Spectroscopy (LIBS)

Laser-induced breakdown spectroscopy (LIBS) offers many attractive analytical advantages with respect to other techniques for the mining industry. LIBS can detect elements from H - Pu, which includes non-metals, such as H, N, F, and O, along with high sensitivity for lighter elements (B, Li, C, K, Ca, Mg, Al, Si, etc.)

### Rapid Composition Monitoring of Raw Li-ion Battery Electrode Materials Using the J200 LIBS Instrument

Laser Induced Breakdown Spectroscopy (LIBS) provides real-time analysis of the raw materials for the Li-ion battery electrode in both powder and pressed pellet form.

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